



CE
0426

MADE IN ITALY



FEMAGINAL PLUS

Ovules

Femaginal Plus ovules is a medical device based on Boric acid, Hyaluronic acid, Polycarbophyl, Lactic acid, Tocopheryl acetate, Vitamin A, 18-Beta-glycyrrhetic acid, Tea tree oil, indicated in the treatment in case of:

 **Vaginal dryness**

 **Irritations**

 **Burning**

 **Infections**

Action of the components:

Lactic acid/Boric acid

IMMEDIATE ACTION

Lowering pH

Boric acid/Lactic acid

Pathogen protection

Polycarbophyl

CONSTANT ACTION


Lowering pH

Hyaluronic acid

Moisturizing and
lubricating action

Vaginal ecosystem: Cyclical variation

Features of the vaginal environment at different ages and physiological status of the woman




Vaginal flora and pH

	PH	LACTOBACILLI	SAPROPHYTIC	PATHOGENS
↑ decrease PH	1	3,8 - 4,5	Minimum Decrease	
	2	4 - 5	Moderate Decrease	Increase
	3	4,5 - 6	Important Decrease	Predominance
FEMAGINAL PLUS	4	> 4,5	Deficit of Production of Lactic Acid	Present
	5	> 5	Reduced amount, unable to produce Lactic Acid	Present
↑ decrease PH	6	> 5,5	Poor, ineffective	Present
	7	6	Absence	Present
	8-10	> 6	Absence	Abundant

Ovules

The importance of the acid pH for the maintenance of a normal vaginal flora



In addition, each ovule contains:

18-Beta-glycyrrhetic acid

Vitamin A and Vitamin E

- Antioxidant activities
- Anti-Decay Cell

Tea tree oil

- Natural preservative
- Without synthetic preservatives

Composition:

One ovule contains: Boric acid, Hyaluronic acid, Polycarbophyl, Lactic acid, Tocopheryl acetate, Vitamin A, 18-Beta-glycyrrhetic acid, Tea tree oil, Semisynthetic triglycerides.

Description: F

Femaginal Plus ovules is a medical device that helps maintain the natural defenses of the vagina and enhances the recovery the normal vaginal flora; containing hyaluronic acid, boric acid and lactic acid it helps delete the symptoms due to unbalanced vaginal hydration and/or vaginal pH. Hyaluronic Acid in its sodium salt form at a 10 mg/ovule dosage exerts an hydrating and lubricating action against vaginal dryness contributes to restore a natural and comfortable vaginal environment. Lactic Acid normally present in the vagina favors an immediate pH correction and facilitates the natural recovery of the delicate vulvo-vaginal ecosystem. Boric Acid, thanks to its persistence in the vaginal environment favors a lowering of pH prolonged over time.

Indication:

Femaginal Plus ovules is useful in the prevention and treatment of vaginal dryness also in mycotic character. It helps reducing irritations, burning and itching.

Dosage:

Insert 1 ovule/day in the vagina – preferably in the evening before sleeping – for 5-7 days after the menstrual cycle or according medical prescriptions.

How to use:

1. Wash vaginal area properly (preferably with both a vaginal wash external and internal);
2. Open the package and after having washed properly your hands, take one ovule from the blister;
3. Gently insert the ovule deep into the vagina, preferably in supine position, avoiding accidental exit of the ovule out of the vagina.

Bibliografy

- 1) Hanna N F, Taylor-Robinson D, Kalodiki-Karamanoli M, Harris J R, McFadyen I R. The relation between vaginal pH and the microbiological status in vaginitis. *Br J Obstet Gynaecol.* 1985;92:1267–1271.
- 2) Valenta C. The use of mucoadhesive polymers in vaginal delivery. *Adv Drug Deliv Rev.* 2005 Nov 3;57(11):1692–712. Epub 2005 Sep 22.
- 3) Delucchi L, Fraga M, Perelmuter K, Cidade E, Zunino P Vaginal lactic acid bacteria in healthy and ill bitches and evaluation of in vitro probiotic activity of selected isolates. *Can Vet J.* 2008 Oct;49(10):991–4
- 4) María Silvina Juárez Tomás, Virginia S. Ocaña, Birgitt Wiese and María E. Nader-Macías Growth and lactic acid production by vaginal *Lactobacillus acidophilus* CRL 1259, and inhibition of uropathogenic *Escherichia coli* *J Med Microbiol* 52 (2003), 1117–1124
- 5) Aroutcheva, A., Gariti, D., Simon, M., Shott, S., Faro, J., Simoes, J. A., Gurguis, A. & Faro, S. (2001). Defense factors of vaginal lactobacilli. *Am J Obstet Gynecol* 185, 375–379
- 6) Boskey, E. R., Cone, R. A., Whaley, K. J. & Moench, T. R. (2001). Origins of vaginal acidity: high D/L lactate ratio is consistent with bacteria being the primary source. *Hum Reprod* 16, 1809–1813
- 7) Ocaña, V. S., Pesce de Ruiz Holgado, A. A. & Nader-Macías, M. E. (1999b). Selection of vaginal H₂O₂-generating *Lactobacillus* species for probiotic use. *Curr Microbiol* 38, 279–284.
- 8) Redondo-López, V., Cook, R. L. & Sobel, J. D. (1990). Emerging role of lactobacilli in the control and maintenance of the vaginal bacterial microflora. *Rev Infect Dis* 12, 856–872
- 9) Melis GB, Ibba MT, Steri B, Kotsonis P, Matta V, Paoletti AM. Role of pH as a regulator of vaginal physiological environment *Minerva Ginecol.* 2000 Apr;52(4):111–21
- 10) Fiorilli A, Molteni B, Milani M. Successful treatment of bacterial vaginosis with a polycarbophil-carbopol acidic vaginal gel: results from a randomised double-blind, placebo-controlled trial. *Eur J Obstet Gynecol Reprod Biol.* 2005 Jun 1;120(2):202–5
- 11) Ceschel GC, Maffei P, Lombardi Borgia S, Ronchi C, Rossi S. Development of a mucoadhesive dosage form for vaginal administration. *Drug Dev Ind Pharm.* 2001 Jul;27(6):541–7
- 12) Maitraie D, Hung CF, Tu HY, Liou YT, Wei BL, Yang SC, Wang JP, Lin CN. Synthesis, anti-inflammatory, and antioxidant activities of 18beta-glycyrrhetic acid derivatives as chemical mediators and xanthine oxidase inhibitors. *Bioorg Med Chem.* 2009 Apr 1;17(7):2785–92. Epub 2009 Feb 21
- 13) Carson C.F. et al. Antimicrobial activity of the major components of the essential oil of *Melaleuca alternifolia*. *J. Appl. Bacteriol.* 78, 264–269, 1995.
- 14) Concha J.M. et al. Antifungal activity of *Melaleuca alternifolia* (tea-tree) oil against various pathogenic organisms. *J. Am. Podiatr. Med. Assoc.* 88, 489–492, 1998
- 15) D'Auria f.d. et al. In vitro activity of tea tree oil against *Candida albicans* mycelial conversion and other pathogenic fungi. *J. Chemother.* 13, 377–384, 2001.

**Thank you for your
attention !**



Health Products with Valuable USP's

www.capricornlifesciences.com.com
info@capricornlifesciences.com.com