


September 17, 2024

Bioanalytical assay strategies and considerations for validation of B cells-depletion assay by flow cytometry

Beata TALAR

*Team Leader, Bioassays Senior Specialist
Mabion S.A.*



AT GLANCE

- > Pharmacodynamics of B-cell-depleting therapeutics is assessed with a combination of direct and indirect methods, including flow cytometry, which is considered a gold standard for measuring B-cell counts in whole blood.
- > Sensitivity is a crucial parameter in bioassays, especially those designed to monitor the pharmacodynamics of cell-depleting therapeutics. Assessing and validating sensitivity involves addressing challenges related to biological variability and lack of biological standards with different levels of dedicated cell populations.
- > "Fit-for-purpose" validation of a flow cytometry assay ensures that the assay is appropriately validated for its intended use, balancing the level of validation effort with the criticality and intended application of the assay.

Introduction

B-cell-depleting monoclonal antibodies and mAb-based products are pivotal in treating various autoimmune and oncological diseases. By selectively targeting and eliminating B-cell populations, these therapies help manage conditions associated with abnormal B-cell function or overactivity, including non-Hodgkin lymphoma, chronic lymphocytic leukaemia, rheumatoid arthritis, systemic lupus erythematosus and multiple sclerosis. Specific depletion with no significant impact on other immune cells is achieved by binding to the CD20 receptor, an antigen expressed primarily on the surface of B-cells. Five drugs from this group are currently marketed across the EU and the US: rituximab (precursor of anti-CD20 mAbs, approved in 1997), ofatumumab, obinutuzumab, ocrelizumab and ublituximab.

Pharmacodynamics of B-cell-depleting therapeutics is assessed via a combination of direct and indirect methods, including flow cytometry and immunohistochemistry monitoring of B-cells, quantitation of serum biomarkers, PK/PD modelling, and, ultimately, evaluation of clinical endpoints. This

multifaceted approach provides a comprehensive understanding of how the tested drug affects target cell populations and overall disease activity, guiding effective treatment strategies and optimizing patient outcomes. Of all listed methods, flow cytometry is considered the gold standard for measuring B-cell depletion and recovery in whole blood after anti-CD20 treatment.

Designing a B-cell depletion assay

One of the most important considerations when designing a B-cell quantification assay is the selection of antigens that will be used to identify the desired population of B lymphocytes. Crucially, the chosen antigen cannot be a therapeutic target, and at the same time, it should be specific to the studied population. Therefore, despite its selective expression, CD20 cannot be used in our assay, as during the anti-CD20 therapy, all receptors are already occupied by administered mAbs, interfering with the binding of fluorescently labelled antibodies.

Like CD20, CD19 receptor is also a B-cell lineage-specific antigen expressed on the cell surface (Figure 1). While CD20 is acquired during the late stages of




Figure 1 Stages of B cells maturations and expression of cell-surface antigens.

B-cell lymphogenesis and is then lost upon differentiation into plasma cells, CD19 expression covers the entire spectrum of early B-cell genesis and maturation.¹ Because CD19 is restricted to B cells, it serves as a specific marker in flow cytometry, allowing for the accurate identification and enumeration of B cells without cross-reactivity with other immune cells.

The B-cell counting assay developed at Mabion utilizes the IVD-certified TBNK reagent (BD Multitest 6-Color TBNK with BD Trucount™ tubes) to determine the percentages and absolute B cell counts in whole peripheral blood. It can also measure the percentages or absolute counts of T lymphocytes (CD3+), natural killer (NK) lymphocytes (CD3-CD16+ and/or CD56+), helper/Inducer T lymphocytes (CD3+CD4+), and suppressor/cytotoxic T lymphocytes (CD3+CD8+).

However, our laboratory customizes the original method to our intended use by:

- ▶ Applying the blood collection tubes that include both anticoagulant (K3EDTA) and cell preservation reagents instead of 'gold standard' K3EDTA tubes prolongs the stability of the cell surface antigens and improves the SSC properties of lymphocytes till the maximum time of 14 days from the time of collection.
- ▶ Testing whole blood subjects derived from RA (Rheumatoid arthritis) and NHL (non-Hodgkin lymphoma) patients in the presence of anti-CD20 treatment.

Following Clinical and Laboratory Standards Institute guideline H62-Validation of Assays Performed by Flow Cytometry, introducing modifications to IVD-certified bioassays requires a systematic and thorough approach to ensure that the modified assay continues to meet performance standards and regulatory requirements.² In this context, the use of blood-stabilizing products is a critical change introduced for the assay. First, cell

preservative compounds can cross-link, masking the target epitopes. Second, these compounds can interfere with tandem fluorochromes, giving false positive results. In this context the By conducting fit-for-purpose validation, our laboratory successfully implemented these modifications, improving assay performance while maintaining results reliability (see section Validation).

TBNK assay performance

The TBNK reagent is composed of seven monoclonal antibodies conjugated to specific fluorochromes (CD45-PrCP-Cy5.5, CD19-APC, CD3-FITC, CD4-PE-Cy7, CD8 APC-Cy7, CD16C-PE D56). The reagent is added to a trucount tube containing a lyophilized pellet of fluorescent beads. Next, the whole blood specimen is directly transferred to the antibodies/beads mix and incubated. During incubation, each monoclonal antibody in the reagent binds to a specific antigen on the surface of blood cells, while beads dissolve, releasing a known number of fluorescent beads. Samples are acquired on a BD flow cytometer (Facs Lyric or FACS Verse) using the predefined protocol on FACS Suite software. The gating strategy of the assay is presented in Figure 2.

After data acquisition, sample gating is reviewed by the analyst and adjusted if required. All predefined protocol changes are automatically recorded in the assay report. Analysts should be guided by Standard Operating Procedure (SOP), which outlines the rules for interpreting and analyzing questionable dot plots. Visual examination of each gate is supported by sample quality questionable error messages, like:

- ▶ Not critical, disease-specific parameters that could also be associated with staining error (%T-sum failure - indicates a large number of double-positive or double-negative T cells, % Lymphosum failure - events could be incorrectly classified as T, B, or NK cells)




Figure 2 A streamlined methodology for PD assessment of anti-CD20 agents: A) TBNK gating strategy*, B) Population hierarchy, C) Identification of blood cells subpopulation on SSC vs. CD45 plot.

- * 1. The software collects 10,000 total events, applies the expert 'lymph' gate on the SSC/CD45 plot, and calculates the number of lymphocytes in the gate. The software collects 10,000 total events, applies the expert 'lymph' gate on the SSC/CD45 plot, and calculates the number of lymphocytes in the gate. The software continues acquiring until the 2500 number of lymphocytes is reached.
- 2. The Truocount beads are gating on SSC/CD19 plot (population all events without lymphocytes – Beads, no Lymph).
- 3. CD3+ lymphocytes are gated on SSC/CD3 plot 5/ CD3+ cells are dividing on CD4+ and CD8+ subpopulation on CD4/CD8 4/ B-cells and NK- defining gates are adjusted on CD19/CD16&56 plot (presenting non-CD3+lymph event).

- ▶ Critical, concerning the acquisition stopping criteria parameters (lymph gate does not contain requested 2,500 events, Insufficient beads acquired (<500).

Quality control samples

Various biological standards are used for flow cytometry, including preserved whole blood, PBMCs (peripheral blood mononuclear cells), and dried leukocytes. However, only a limited number of standards are provided with certified target values, such as absolute cell counts. If they are certified, the range for target value is extensive, and often, it is not

provided for each subset analyzed with a given method. Consequently, they cannot be considered 'real' reference standards. The BD™ Multi-Check Control and CD-Chex Plus (Streck) are stabilized human leukocytes and erythrocytes in a preservative medium dedicated to complete process control for TBNK immunophenotyping of whole blood samples. There is a control for antibody staining, red blood cell (RBC) lysis, instrument setup and performance, and data analysis on FACS. The experimental cell number value of the specific lymphocyte subpopulations must be within limits specified by the manufacturer in the certificate of analysis.

Following the manufacturer's recommendations, the biological standard should be analyzed at least once per analytical day. However, some laboratories, including ours, analyze the biological standard parallel to the test samples during each analytical run.

Day-to-day instrument and assay standardization, monitoring and setup

PQC (Performance Quality Control) with CS&T beads (Cytometer Setup and Tracking beads) are used to define the baseline performance of the cytometer. They check laser alignment, measure %rCV, linearity, resolution and laser power, compare PMT voltages to characterization QC (baseline), and update spillover values and PMTVs for default settings. Daily measurements are automatically entered into Levey-Jennings to detect increased random errors and shifts or trends in cytometer performance. This meticulous monitoring of FACS performance measurements over time ensures consistent instrument performance, promptly addresses issues and maintains the quality and reliability of data.

Validation

The CLSI H62 document marks a significant advancement in standardizing flow cytometry assay

validation practices, building on the foundational consensus papers published over the last decade.^{7,8} Before its publication, no official guidelines addressed the validation requirements of flow cytometric methods used in drug development or clinical testing.

Due to the complexity of flow cytometry technology, lack of certified reference materials, and different categories of bioanalytical data (qualitative, quantitative, quasi-quantitative), not every validation parameter can be assessed.⁸ Assay specificity, sensitivity (limits of detection /quantification), precision (repeatability), stability and reference ranges are always validated for the flow cytometry method. Linearity, standard calibration, and interference (matrix, drug) can also be validated in some cases. However, accuracy, selectivity, range of quantification, incurred sample reanalysis, normal signal distribution, and prozone effect usually cannot be evaluated in a direct manner. Furthermore, the specific validation parameters to be assessed depend on the intended data use. If it changes during the assay's lifetime, additional validation parameters can be addressed at that time. This iterative approach is known as "Fit-for-Purpose" validation.²

The TBNK reagent is used for the measurement of [%] and absolute counts of Ts, Th, B, and NK Cell.

Table 1 Quality Control of Instrument Performance and Instrument Setup

Activities	Frequency	Reagents	Description
PQC	Daily	CST beads	Performance quality control and instrument set up
TBNK ASSAY SET UP	Daily	CST beads	Assay set up
CQC	Every 6 months	CST beads	Characterization QC that, along with additional more detailed analyses and an automatic laser alignment
REFERENCE SETTING RESET AND ACTUALIZATION	Every 60 days	FC beads	To measure and reset instrument spectral overlap values

Table 2 Assessment of validation parameters for quasi-quantitative determination of B-lymphocytes in human whole blood in the presence of anti-CD20 treatment

Parameter	Our analytical approach
MATRICES	Disease-state samples or initial validation with samples from healthy individuals, followed by confirmation of assay precision, sensitivity, and selectivity using disease-state samples.
PRECISION	A multifactorial design is developed to evaluate multiple factors in a single setting, allowing for intra-assay precision, inter-assay precision, inter-operator precision, inter-instrument precision, and reagent lot-to-lot consistency.
ACCURACY	Data obtained from flow cytometry are considered quasi-quantitative because the numerical result generated from the test sample (cell/ μ L) is proportional to the sample (truecount tubes with a certified number of beads) but not generated from a calibrator or reference material, so true accuracy was not determined during the method validation. Furthermore, we could not determine if our method is accurate because there is no reference standard containing a certified number of B cells (at different cell count levels). Although target ranges of BD Multicheck Control and CD Check (Streck) are provided, the ranges are extremely broad and thus of little value for accuracy assessment. As an alternative to true accuracy measurements, results from the inter- and intra-assay precision experiments should be compared with the reference laboratory after instrument-to-instrument standardization.
SELECTIVITY/ MATRIX EFFECT	<ul style="list-style-type: none"> > Individual highly lipemic, hemolyzed, clotted and disease-state samples are used to verify the appropriate subpopulation separation. > Evaluation % of the difference for B absolute counts results obtained for subjects collected to blood collection tubes that include both anticoagulant (K3EDTA) and cell preservations vs. 'gold standard' K3EDTA tubes.
SENSITIVITY	For functional LoQ determination, low-B-cell-concentrated samples are created by diluting whole blood using autologous plasma reconstituted with RBCs to achieve CD3-/CD19+ absolute counts of 5, 10, 20, 40 (cells / μ L) (LoQ candidates). These absolute B cell count values must be inside the analytical measurement range specified by the manufacturer (5-2000 cells / μ L). Three independent repetitions (three separate runs) should be prepared from each of the concentration pools of 5 test samples.
STABILITY	<ul style="list-style-type: none"> > Post-collection stability - age of blood (AOB) refers to the time between the blood draw and the time of mixing the blood with the reagents. > Post-sample processing stability: <ul style="list-style-type: none"> - Post-staining stability - Age of stain (AOS) refers to the incubation time of blood with antibodies in Trucount tubes. - Post-lysing stability - Age of fixation (AOF) refers to the time between the end of the lysing step and the start of specimen acquisition. > Thermal cycles stability
CARRYOVER	Carryover is evaluated as % Lymphocyte carryover from Stained Sample to Blank sample should be less than or equal to 0.5% when approximately 2500 Lymph events are acquired.

However, the CD3-CD19+ B cells are enumerated for pharmacodynamics evaluation of anti-CD20


therapeutics, and consequently, the results obtained for this subpopulation are considered critical for

determining assay validity. This assay is classified as quasi-quantitative because the absolute cell count can be determined by comparing cellular events to the certified number of bead events. In Table 2 on the next page, we present our approach to validating the quasi-quantitative determination of B lymphocytes in human whole blood in the presence of anti-CD20 treatment. For each validation parameter, the acceptance criteria were predefined before the start of validation.

Since we modified the original test protocol using blood-stabilizing products, monitoring the correct separation of the TBNK subpopulation over time was critical during the validation phase. Ageing of test

samples can cause both losses of TBNK antigen detection (due to protein degradation) and change in Leucocytes' morphology (loss of the difference in SSC properties). Both cases could affect the absolute lymphocyte counts' measurements and inclusion of basophils (less bright CD45, low side scatter) or monocytes (less bright CD45, moderate side scatter) in a lymphocyte cluster.

Our validation data showed that test samples collected in blood collection tubes that include both anticoagulant (K3EDTA) and cell preservation reagents (P2) remain stable at 2-8°C up to 14 days from the day of collection and at RT up to 3 days from the time of collection compared to the baseline



Abbreviations: P1 – gold-standard K3EDTA tubes without preservatives; P2 - blood collection tubes that include both anticoagulant (K3EDTA) and cell preservation reagents, Time of testing: T0 ($\leq 6h$ from the draw), T24 ($24h \pm 2h$ from the draw), T7d ($7days \pm 1day$ from the draw), T10d ($10 days \pm 1day$ from the draw), T14d ($14 days \pm 1day$ from the draw).

Figure 3 Cell preservatives effect on TBNK subpopulations separation.

The visualization of the lymphocytes tested via cell markers CD45, CD3, CD19, CD16+56 for the control (P1 and P2 < 6 hours, Day 0), P2 after 24h, on day 7, and 14, have shown that the Lymphocytes from P1 samples exhibited higher side scatter (SSC) properties than the equivalent cells in the P2 tubes; however, the SSC properties also increase in time (from T0 to 14days) in P2 tubes. On the other hand, the fluorescence characteristics are mostly stable in both time and storage conditions. These results demonstrate that although stabilization affects the cells' physical properties, the different cell populations can be adequately identified.

sample (collected in gold-standard K3EDTA tubes without preservatives (P1)), as presented in Figure 3. Furthermore, leukocyte subpopulations are adequately separated at any tested point in time.

Figure 3 shows the lymphocytes' visualization of DonorX tested via cell markers CD45, CD3, CD19, CD16+56 for the control (P1 and P2 tested <6 hours after drawing, Day 0), P2 after 24h, on day 7, and 14. Lymphocytes from P1 samples exhibited higher side scatter (SSC) properties than the equivalent cells in the P2 tubes, and the SSC properties also increased in time (T0 to 14 days) in P2 tubes. However, the leucocyte subpopulations are adequately separated at any point in time.

After validating an assay with a comprehensive QC process, it's crucial to continuously monitor and maintain performance to ensure it aligns with the established specifications.

TBNK immunophenotyping assay in the monitoring and management of oncological therapies

TBNK immunophenotyping is an invaluable tool in managing oncological and autoimmune disorders. Developed initially to monitor CD4 T-cell counts in HIV-positive patients⁹, this assay has evolved significantly, finding critical applications across various medical fields. The ability to quantify and characterize T-cell, B-cell, and NK cell populations provides profound insights into immune system dynamics, aiding in the diagnosis, monitoring, and therapeutic assessment of numerous diseases.

TBNK immunophenotyping can be used for the following purposes:

- ▶ **Real-time Monitoring:** Provides immediate feedback on the immune system's response to therapy, allowing for rapid adjustments in treatment strategy.
- ▶ **Mechanistic Insights:** Helps in understanding how drugs modulate the immune system, providing insights into their mechanisms of action.
- ▶ **Predictive Value:** Early changes in TBNK profiles can predict long-term outcomes, helping to identify responders and non-responders.
- ▶ **Safety Management:** Identifies potential immune-related adverse events early, allowing for proactive management of toxicities.

Specific examples of TBNK assays' application include:

- ▶ **Checkpoint Inhibitors (e.g., Anti-PD-1/PD-L1, Anti-CTLA-4)**¹⁰
 - Pembrolizumab (Keytruda) and Nivolumab (Opdivo): These anti-PD-1 antibodies are used in the treatment of melanoma, renal cell carcinoma, and non-small cell lung cancer. TBNK immunophenotyping can track the changes in T-cell populations, significantly increasing the number of effector T-cells and decreasing the number of regulatory T-cells, which correlate with therapeutic response.
 - Ipilimumab (Yervoy): An anti-CTLA-4 antibody used in the treatment of melanoma and other solid tumors. Monitoring the increase in activated T-cells and reduction in regulatory T-cells helps assess treatment efficacy and manage immune-related adverse events.
- ▶ **CAR-T Cell Therapy (e.g., Tisagenlecleucel, Axicabtagene Ciloleucel)**
 - Tisagenlecleucel (Kymriah): Indicated for B-cell acute lymphoblastic leukaemia (ALL) and diffuse large B-cell lymphoma (DLBCL). TBNK immunophenotyping is crucial to track the expansion and persistence of CAR-T cells in the patient's blood, as well as the depletion of B-cells, which indicates the CAR-T cells' effectiveness.


- Axicabtagene Ciloleuceel (Yescarta): Used for treating DLBCL. Monitoring involves tracking CAR-T cell populations and their impact on normal B-cell counts to manage potential B-cell aplasia.
- ▶ Bi-specific T-cell Engagers (BiTEs) (e.g., Blinatumomab)
 - Blinatumomab (Blinicyto): Used in the treatment of B-cell precursor ALL. TBNK immunophenotyping helps in monitoring the engagement between T-cells and B-cells, indicating the drug's efficacy in redirecting T-cells to kill B-cell leukaemia cells.
- ▶ Monoclonal Antibodies (e.g., Rituximab, Daratumumab)
 - Rituximab targets CD20 on B-cells and is used as a treatment for non-Hodgkin lymphoma (NHL) and chronic lymphocytic leukaemia (CLL). TBNK immunophenotyping assesses B-cell depletion and recovery post-therapy.
 - Daratumumab (Darzalex): Targets CD38 on plasma cells and is used for treating multiple myeloma. Monitoring involves tracking the depletion of CD38-positive cells and changes in other lymphocyte populations.
- ▶ Immune Modulators (e.g., Lenalidomide)
 - Lenalidomide (Revlimid): Used as a treatment for multiple myeloma and mantle cell

lymphoma. TBNK immunophenotyping helps evaluate the drug's effect on enhancing T-cell and NK cell activity while reducing regulatory T-cells, contributing to its anti-tumour effects.

Summary

Flow cytometry assays have a variety of intended purposes that include basic research, drug discovery, and the exploratory, secondary, and primary endpoint of clinical trials. In designing a fit-for-purpose validation for a B cell-depleting assay, it is crucial to tailor the validation experiments to meet the specific requirements of the assay's intended use. This approach ensures that the assay delivers reliable, "accurate" results and is efficient and practical. Several parameters, such as accuracy, sensitivity, and selectivity, pose significant challenges in validating a TBNK assay, especially when dealing with quasi-quantitative data. This paper describes the solutions employed in our laboratory to address these challenges effectively. Remember:

- ▶ Any modifications made by the laboratory to the assay manufactured as IVD's instructions will render the test no longer FDA-cleared/approved. The laboratory must verify/validate these modifications in terms of the intended use.
- ▶ In clinical trials, the use of solutions that effectively extend the 'lifetime' of the test sample is not just a necessity, but a critical

 provides **comprehensive flow cytometry services** tailored to the needs of biopharmaceutical companies, clinical laboratories and research institutions. Our experienced scientists offer robust support for various flow cytometry endpoints using various cell types, including primary cells and human or animal-origin cell lines. Whether you require non-regulated or regulated assays (GMP/GLP/GCP), we can provide phase-appropriate flow cytometric bioanalysis and confidently assist you utilizing the following platforms: BD FACSLyric™ and BD FACS Verse™ powered by the BD FACSuite™ Acquisition and Analysis Application.

need, given the time-sensitive nature of the trials and the need for reliable data.

- ▶ Validation in flow cytometry is not straightforward; each validation parameter's experiment is designed based on the purpose of the test. There is no universal design for all methods, but we have developed a specific approach that we have found to be effective. You can refer to Table 3 for a detailed description of our specific approach.

If you are interested in this publication or would like to learn more about flow cytometry methods, please get in touch with us to collaborate.

References

1. Edwards JC, C. G. (2006, May). *B-cell targeting in rheumatoid arthritis and other autoimmune diseases*. *Nat Rev Immunol*.
2. CLSI, Clinical and Laboratory Standards Institute. (2021). *H62 Validation of Assays Performed by Flow Cytometry*.
3. der Strate BV, L. R. (2017). *Best practices in performing flow cytometry in a regulated environment: feedback from experience within the European Bioanalysis Forum*. *Bioanalysis*.
4. Green CL, B. L. (2011). *Recommendations for the validation of flow cytometric testing during drug development: I instrumentation*. *Journal of immunological methods*.
5. Hara DM, X. Y. (2011). *Recommendations for the validation of flow cytometric testing during drug development: II assays*. *Journal of immunological methods*.
6. Johansson U, B. D. (2014). *Guidelines on the use of multicolour flow cytometry in the diagnosis of haematological neoplasms*. *British Committee for Standards in Haematology*. *British journal of haematology*.
7. L, S. (2018). *2018 White Paper on Recent Issues in Bioanalysis: focus on flow cytometry, gene therapy, cut points and key clarifications on BAV (Part 3 - LBA/cell-based assays: immunogenicity, biomarkers and PK assays)*. *Bioanalysis*.
8. Selliah N, N. V. (2018). *Flow Cytometry Method Validation Protocols*. *Curr Protoc*.
9. Omana-Zapata I, Mutschmann C, Schmitz J, Gibson S, Judge K, Aruda Indig M, Lu B, Taufman D, Sanfilippo AM, Shallenberger W, Graminske S, McLean R, Hsen RI, d'Empaire N, Dean K, O'Gorman M. *Accurate and reproducible enumeration of T-, B-, and NK lymphocytes using the BD FACSLyric 10-color system: A multisite clinical evaluation*. *PLoS One*. 2019 Jan 28;14(1):e0211207.
10. Carlino MS, Larkin J, Long GV. *Immune checkpoint inhibitors in melanoma*. *Lancet*. 2021 Sep 11;398(10304):1002-1014. doi: 10.1016/S0140-6736(21)01206-X. PMID: 34509219.