



Maggie®

Two Step Cartridge-Based Solution



Volume 

Viscosity 

Device Type

Primary Container

Dosage Type

Route of Administration

Activation Method

Injection Feedback

Disposable autoinjector

3.0 mL standard ISO cartridge with pre-installed needle (NIT®)

Fixed Dose | Optional automated priming and titration

Subcutaneous | Intramuscular

Needle cover

Continuous audible clicks during injection, viewing window

Handling Instructions



1 Hold device with cap upright. Unscrew cap until needle cover automatically extends. Then remove cap.



2 Press and hold against injection site to activate.



+ Discard in sharps container.

Key Features

- ✓ Simple 2-step operation
- ✓ Needle cover activation
- ✓ Continuous audible and visual feedback
- ✓ One-handed operation
- ✓ Permanently hidden needle
- ✓ Compact design
- ✓ Customizable needle length and gauge
- ✓ Passive locking needle cover
- ✓ Customized to client specifications



SHL MEDICAL

The SHL Medical logo and Maggie are registered and/or non-registered trademarks of SHL Medical AG or its affiliates, in the United States, Europe, and other jurisdictions. Final product specifications are subject to change based on customer requirements. Contact info@shl-medical.com for more information.

www.shl-medical.com