

Volume Viscosity 

Device Type

Primary Container

Dosage Type

Route of Administration

Activation Method

Injection Feedback

Disposable autoinjector

5.0 mL cartridge with pre-installed needle (NIT®)

Fixed dose | Optional automated priming and titration

Subcutaneous | Intramuscular

Needle cover

Continuous audible clicks during injection, viewing window

## Handling Instructions



Hold device with cap upright. Unscrew cap until needle cover automatically extends. Then remove cap.



Press and hold against injection site to activate.



Discard in sharps container.

## Key Features

- ✓ Simple 2-step operation
- ✓ Needle cover activation
- ✓ Continuous audible and visual feedback
- ✓ One-handed operation
- ✓ Permanently hidden needle
- ✓ Compact design
- ✓ Customizable needle length and gauge
- ✓ Passive locking needle cover
- ✓ Customized to client specifications



The SHL Medical logo and Maggie are registered and/or non-registered trademarks of SHL Medical AG or its affiliates, in the United States, Europe, and other jurisdictions. Final product specifications are subject to change based on customer requirements. Contact [info@shl-medical.com](mailto:info@shl-medical.com) for more information.

[www.shl-medical.com](http://www.shl-medical.com)