



i+Med 
BIOMEDICAL
INTELLIGENCE

BtH Dermal
Fillers

BtH Dermal Fillers

It is a line of hyaluronic acid-based dermal fillers specially designed to smooth wrinkles, plump lips, restore skin volume and redefine the facial contour.

The gels are manufactured with hyaluronic acid of non-animal origin, obtained by biofermentation and highly purified. Each product of the line has been developed bearing in mind a precise treatment objective.

WHAT IS BtHLINE

BtHLINE is intended as a dermatological treatment for aesthetic purposes. It is formulated with linear HA of high molecular weights.

Administered by intradermal injection, the product is indicated to complement the extracellular matrix and intradermal tissue, restoring the anatomical structures of the skin and promoting its hydration, revitalization and firmness.



BtHLINE

BtHLINE is presented in a pre-filled syringe of 1mL, ready for use. It is supplied in a single blister together with two sterile needles (30G) for single use, as accessories.

MODEL	BtHLINE 1,4%
VISCOSITY (cP)*	10.000 - 40.000
MOLECULAR WEIGHT	High Molecular Weight.
COMPOSITION	Hyaluronic acid (non-animal origin)
ENDOTOXIN CONTENT	<0,5 EU/mL
SIZE	1mL

* Viscosity measured at shear rate (SR) 1s at 25°C. Measured after sterilizing.

BtHLINE PROPERTIES


- + Sterile product
- + Biocompatible product
- + Thanks to its exclusive compositions it allows you to achieve a natural look
- + Restores aged skin structures.
- + Improves skin hydration.
- + Revitalizes skin and improves its firmness.



WHAT IS BtHCROSS

BtHCROSS is a biocompatible hydrogel of cross-linked hyaluronic acid obtained by biofermentation and highly purified.

BtHCROSS has been specifically design for dermatological treatment with aesthetic purposes. Administered by intradermal injection, the product is indicated to restore and correct facial volume and fill in wrinkles and folds.



BtHCROSS

BtHCROSS is available in three models:

i+Med has developed three different models of cross-linked dermal fillers to offer the most adapted solutions depending on the needs of each patient.

MODEL	BtHCROSS 1,5%	BtHCROSS 1,75%	BtHCROSS 2%
VISCOSITY (cP)*	20.000 - 30.000	30.000 - 40.000	50.000 - 60.000
MOLECULAR WEIGHT	High Molecular Weight.		
COMPOSITION	Cross-linked Hyaluronic acid (non-animal origin)		
ENDOTOXIN CONTENT	<0,5 EU/mL		
SIZE	1mL		

* Viscosity measured at shear rate (SR) 1s at 25°C. Measured after sterilizing.

BtHCROSS PROPERTIES

- + Sterile product.
- + Biocompatible product.
- + Thanks to its exclusive compositions it allows you to achieve a natural look.

Depending on the area to be treated and the need of the patient, there are 3 different models that provide the best solution of each need:

BtHCROSS 1,5%	<ul style="list-style-type: none">• Corrects superficial wrinkles such as perioral lines.• Enhances and profiles lips.
BtHCROSS 1,75%	<ul style="list-style-type: none">• Corrects medium to deep skin depressions such as nasolabial folds.
BtHCROSS 2%	<ul style="list-style-type: none">• Remodels and restores volume of the facial contour.



PARQUE TECNOLÓGICO
DE ÁLAVA

Calle de Albert Einstein, 15
01510 Vitoria-Gasteiz
(Álava-Spain)

(+34) 945 56 11 34
imasmed@imasmed.com

www.imasmed.com