

About Us

When it comes to the pharmaceutical industry, trust and quality are paramount. Since our inception in 2013, Gelcaps Industries has proudly stood as a beacon of reliability, consistently delivering the highest quality empty hard gelatin capsules to meet the stringent demands of the healthcare sector.

PRODUCTS

TSE/BSE-Free Capsules

We make EDQM (European Directorate for the Quality of Medicines & Healthcare) certified and CEP grade hard gelatin capsules made from Type B Bovine Bone. These special variants are made in order to prevent the risk of Transmissible Spongiform Encephalopathy (TSE) and Bovine Spongiform Encephalopathy (BSE).

XPress-Fit Tabsules

Unlike regular tablets, these tables are made from pressing flexible caps over the tablets using a cold-shrink process. They are more attractive, uncommon, odorless, easier to swallow and faster to digest. These XPress-Fit tabsules can give your product a new unique look without changing any formulation or color.

Pearl Capsules

The capsules made from a typical silver shiny pigment are called 'pearl capsules'. Their attractive color shades make them unique and act as a product differentiator. These types of capsules are made on special demand from certain companies based on their marketing needs.

Preservative-Free Capsules

The empty hard gelatin capsules manufactured at our plant are made without using any preservatives making them safer to consume. These variants are more natural and easier to dissolve in the patients' body.

SLS Free Capsules

These capsules are specially manufactured without using SLS (Sodium Lauryl Sulphate) into manufacturing of empty hard gelatin capsules.

Halal/Kosher Certified Raw Material

Respecting the religious laws and beliefs of every religion, we make our gelatin capsules only from selective cattle bones that are permitted under Halal/ Kosher laws.



The hallmarks of our quality



Product Details

- 100% compliance to regulatory market
- WHO-GMP certified
- Microbiological Precision: Our software-monitored incubators perform rigorous microbiological quality checks in accordance with harmonized monographs.
- Pristine Environment: Active air sampling units maintain an immaculate environment at every production stage.
- Stability Guaranteed: We adhere to ICH guidelines with dedicated stability chambers, meticulously monitored by software.
- Light Sensitivity Testing: Our chambers also facilitate photostability testing for light-sensitive colorants.
- Precision Checks: We conduct thorough assessments of dimensions, attributes, and suitability of products during the filling machine operation.
- Ultimate Protection: To safeguard against transit damage, our final products are meticulously packed in 110-micron antistatic food-grade polythene bags before being batch packed in 5-ply pre-printed corrugated boxes.



SUMMARY OF DIMENSIONS FOR EMPTY HARD GELATIN CAPSULE SHELLS

CAPSULE SIZE	OEL	0	1	2	3	4
1.Capacity						
Body Volume (ml)	0.78	0.68	0.5	0.37	0.3	0.21
2. Weight						
Targeted Wt In (mg)	104	96	76	63	50	39
Wt. Tolerance In (%)	Tolerance On Avg is ± 5% & Individual Capsule ± 10%					
3.Empty Overall Length						
Cap Nominal Length (mm)	11.7	10.7	9.8	9	8.1	7.2
Body Nominal Length (mm)	20.3	18.5	16.6	15.2	13.6	12.2
Nominal Closed Joined Length (mm)	23.2	21.4	19.3	17.8	15.8	14.4
Length Tolerance In (mm)	0.4	0.4	0.4	0.4	0.4	0.4
4. Outside Diameter (Caps & Body)						
Cap Nominal Diameter (mm)	7.66	7.66	6.93	6.38	5.85	5.33
Body Nominal Diameter (mm)	7.34	7.34	6.63	6.09	5.57	5.06
Diameter Tolerance In (mm)	0.03	0.03	0.03	0.03	0.03	0.03

The hallmarks of our quality



GELCAPS INDUSTRIES

ADDRESS

Corporate Office: 21st Floor, Galaxy Tower, Hyderabad Knowledge City, TSIC Raidurgam Panmaktha, Rangareddy District, Hyderabad-500081, Telangana, INDIA

FACTORY:

Plot No. 11, Block-B, APIIC Industrial Park, Menakur (V), Naidupet (M), SPSR Nellore Dist- 524126, Andhra Pradesh, INDIA.
 ☎ Ph: 07731950022