SAFETY ⇔ CLEANROOM PRODUCTS

GREEN GHD STOPPER BOWL COVERS



Save time by using the stopper bowl cover instead of wrapping with sheets of material and autoclave tape. When used in conjunction with a secondary autoclave bag, it creates an effective wrapping system. After the outer bag is removed, this stopper bowl cover should remain in place, providing protection from microbial and particulate contamination while the remaining set-up activities are performed.

TYPICAL APPLICATION & CHARACTERISTICS:

- Steam Sterilization Indicator
- Autoclavable up to 135°C
- Incremental improvement compared to cutting and wrapping with Bioshield™ and autoclave tape
- Eliminates need to cut material and the use of particulate generating autoclave tape
- Fast and consistent
- Creates secondary level of protection during aseptic presentation and setup of equipment
- Avoid sterilization Re-validation
- Manufactured and packaged in ISO8 certified Cleanroom
- Conveniently packaged: 1 per heat-sealed labeled poly bag
- Certificates of Conformance sent with all shipments
- Made in USA

MATERIAL SPECIFICATIONS:

PROPERTY	UNITS	GHD MATERIAL
Basis Weight	Gm/m2	57
Dry Tensile Strength (MD)	grams	8,100
Dry Tensile Strength (CD)	grams	7,100
Elongation Dry	%	27.3
Air Permeability	liters / 100 sq. cm	110
Suter Water Head (pressure test)	mm	400
Dry Mullen Burst	Grams / sq.cm	2,150
Flammability	seconds	8

GREEN GHD STOPPER BOWL COVERS



KEYSTONE PART NUMBER	DIMENSIONS	
020GHD	Fits 27.94 cm to 40.64 cm Diameter	
024GHD	Fits 30.48 cm to 50.8 cm Diameter	
030GHD	Fits 50.8 cm to 66.04 cm Diameter	
036GHD	Fits 66.04 cm to 76.2 cm Diameter	
042GHD	Fits 71.12 cm to 88.9 cmDiameter	
048GHD	Fits 76.2 cm to 101.6 cm Diameter	
060GHD	Fits 101.6 cm to 127 cm Diameter	

